**TROUSSE DE MONITORAGE DE PRESSION
SOUS-DURALE POST-CRANIOTOMIE
MODELE 110-4G**

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**POST CRANIOTOMY
SUBDURAL PRESSURE
MONITORING KIT**

MODEL 110-4G

DIRECTIONS FOR USE

Read carefully prior to use

WARNING:

This product contains a chemical (Ethylene Oxide) known to the State of California to cause cancer, birth defects or other reproductive harm.

FEDERAL (U.S.A.) LAW RESTRICTS THIS DEVICE TO SALE
BY OR ON THE ORDER OF A PHYSICIAN.


 **INTEGRA™**
Neuro Sciences
CAMINO®

*POST CRANIOTOMY SUBDURAL PRESSURE MONITORING KIT

MODEL 110-4G

*The Post Craniotomy Subdural Pressure Monitoring Kit was developed in cooperation with Richard C. Ostrup, M.D., Thomas Luerssen, M.D., and Lawrence F. Marshall, M.D. (San Diego).

Figure 1

SYSTEM DESCRIPTION

The Post Craniotomy Subdural Pressure Monitoring Kit contains a sterile transducer-tipped pressure monitoring catheter to be used as a diagnostic tool for rapidly determining and continuously monitoring intracranial pressure. Since its method of measuring pressure is unique, please read this section carefully.

In addition to the catheter, the following accessories are included with the 110-4G (see Figure 1).

- Zero adjustment tool.
- (1) each 0.125" (3.17 mm) stainless steel trocar with scalp tunneling sheath.

Unlike ordinary pressure monitoring systems, the Camino® Catheter has a miniature transducer at the distal tip. This unique design eliminates the need for a "fluid-filled system" to carry pressure waves to an external transducer. Thus, the problems inherent in such fluid-filled systems are eliminated. The result is a precise pressure measurement and an artifact free, high fidelity waveform trace.

Note: A hand drill and various standard surgical instruments and supplies are required to place the ICP monitoring system.

A complete set of instruments and supplies is available from Integra NeuroSciences in our Cranial Access Kits.

INDICATIONS

The use of the Post Craniotomy Subdural

Pressure Monitoring Catheter by a qualified neurosurgeon is indicated when direct pressure measurement in the subdural space, post craniotomy, is clinically important.

CONTRAINDICATIONS

This device is not intended for any use other than that indicated.

RISK AND COMPLICATIONS

Maintenance of sterility during placement and subsequent handling of the Post Craniotomy Subdural Pressure Monitoring Kit is essential. Sterile technique should be used at all times when inserting, adjusting and securing the Camino Catheter.

If monitoring is continued for more than 5 days, placement of a new system under sterile conditions is recommended.

Placement of the Post Craniotomy Subdural Pressure Monitoring Catheter must be carried out by a qualified neurosurgeon using standard surgical procedure and skill. Appropriate measures to avoid infection and complications are the sole responsibility of the neurosurgeon in charge.

PRECAUTIONS

- Extreme bending and/or kinks can impair the performance of the Fiber Optic Pressure Transducer. Exercise caution when handling the Catheter.
- The Catheter is designed for **SINGLE USE ONLY. DO NOT RESTERILIZE OR REUSE.**

Camino 110-4G is supplied in a single use package and is guaranteed to be sterile and non-pyrogenic unless opened or damaged. The product is intended for use as an intracranial pressure monitor and is not to be reused. Any attempts to resterilize or reuse the product will damage the catheter and impair its ability to function as intended.

- Use aseptic technique throughout procedures.
- Maintain the insertion site with regular, meticulous redressing using aseptic technique.
- Do not attach anything to transducer air vent. Vent must remain open for proper operation (Figure 1).

CAUTION

Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

INSTRUCTIONS FOR USE

The Camino Catheter can be used with the Camino 420/V420 series or with Camino MPM series pressure monitors. For Camino monitor set-up and use, refer to the appropriate Camino monitor Directions For Use.

CATHETER PREPARATION PRIOR TO INSERTION

Remove the Camino Catheter from its sterile package and firmly attach the transducer connector to the preamp connector (Figure 2). If the Camino display does not read zero after a short system self check display, use the tool from the catheter kit to turn the zero adjustment on the bottom side of the transducer connector until the Camino display reads zero (Figure 3).

To measure post-craniotomy subdural pressure, first prepare the Camino transducer-tipped catheter as previously instructed. Choose the burrhole through which the Camino catheter will enter and notch the corner. After the catheter has been zeroed, position the tip of the catheter within the trocar's guide tube, and use the trocar to tunnel the catheter under the scalp, toward the craniotomy site. Then remove the trocar.

Place the tip of the catheter on the brain tissue, under the dura, opposite the notched burrhole. Locate the tip under the intact skull, to isolate it from local pressures which may be caused by movement of the bone flap. (Figure 4).

Close and suture the dura using standard neurosurgical procedure. Position the catheter through the notched burrhole, and replace the bone flap. Suture around, not through, the catheter, and secure it to the scalp.

After the site has been dressed, secure the preamp connector to the patient to protect the fiberoptic catheter.

IT IS RECOMMENDED THAT THE CATHETER BE DISCONNECTED FROM THE PRE-AMP CONNECTOR WHEN THE PATIENT IS MOVED. THIS WILL NOT AFFECT CALIBRATION.

CONTINUOUS PRESSURE MONITORING

Since the Camino Catheter has a miniaturized transducer at the distal tip, it requires no fluid-filled system. Thus, the need for an external transducer, pressure tubing is eliminated. As a result, pressure may be monitored continuously without recalibration.

Figure 2


Figure 3



Figure 4

INTERNATIONAL MARKINGS

CATHETER EXTENSION CABLE CONNECTIONS

Symbol also signifies:

TYPE BF EQUIPMENT: Protected against electric shock, having an F-type isolated (floating) applied part.

Sterilized using ethylene oxide.

Date of expiration.

Do not reuse, single use device.

Manufacturer.

Do not resterilize.

TYPE CF EQUIPMENT: Provided with a degree of protection against electric shock higher than Type BF equipment, having an F-type isolated (floating) applied part; intended for direct cardiac application.

Refer to Directions for Use.

Lot Number.

Rx ONLY

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

Authorized representative in the European Community.

SPECIFICATIONS

Transducer size	4FR (1.35mm Diameter)
Transducer type	Fiber optic
Frequency response (system)	
Model MPM-1/SPM-1	100Hz (-3dB)
Model V420/M420/420XP	120Hz (-3dB)
Model 420	33Hz (-3dB)
Measurement range (system)	-10 to 125 mmHg
Zero drift (system)	
First 24 hours (maximum)	0 ± 2 mmHg
5 days (typical)	less than ± 1 mmHg per day
Temperature coefficient	Max of 3 mmHg over temperature range of 22°C to 38°C (70°F-100°F)
Linearity and hysteresis (system)	
Pressure range:	
-10 to 50 mmHg	± 2 mmHg or better
51 to 125 mmHg	± 6% of reading or better
Reference pressure	Atmosphere
Overpressure	-700 to 1250 mmHg

Camino and Integra NeuroSciences are registered trademarks at Integra LifeSciences Corporation in the United States and/or other countries. The Integra wave logo is a trademark of Integra LifeSciences Corporation © 2010 Integra LifeSciences Corporation. All Rights Reserved.

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Integra NeuroSciences Ltd.,